

Centrifuge protocol & PRP preparation

సెంట్రీఫ్యూజ్ & PRP ప్రిపరేషన్



TODAY'S OBJECTIVES · ఈ రోజు లక్ష్యాలు

- ✦ Operate the centrifuge correctly and safely
- ✦ Prepare PRP to the clinic's exact protocol
- ✦ Maintain sterility from tube to syringe

THEORY · థియరీ

Centrifuge basics

Always balance: tubes of equal weight opposite each other; add a water-filled balance tube if you have an odd number. Close the lid, set speed and time exactly as protocol, never open before it stops completely. Typical single-spin hair PRP: around 1500 rpm for 8–10 minutes, or the brand's stated protocol — never improvise.

Separating the PRP

After the spin: plasma on top, buffy coat line, red cells at the bottom. Using a sterile syringe and long needle or the kit's transfer device, draw the plasma just above the buffy coat, including the buffy layer as per protocol, without disturbing the red cells. Aim for the volume the doctor asked for; typical scalp session 4–6 ml.

Sterility chain

Sanitise and glove, work on a sterile draped tray, open kits without touching inner surfaces, never leave a prepared syringe uncovered, use it within the protocol time window, and label with the client's name. Discard immediately if there is any doubt about contamination — a scalp infection is far more expensive than a kit.

👋 Hands-on today

- ✦ Balance and run the centrifuge 3 times with dummy tubes
- ✦ Prepare PRP under supervision from a real draw and have the volume checked
- ✦ Set up and break down a sterile tray, audited

⚠️ Safety & contraindications

- ✦ Never open a running centrifuge; never run unbalanced
- ✦ One client = one kit; never share or reuse anything

KEY POINTS TO REMEMBER · గుర్తుంచుకోవలసినవి

- ✦ Balance every spin
- ✦ Follow brand rpm and time exactly
- ✦ Typical scalp dose 4–6 ml
- ✦ Label and use within the time window

📖 Homework

Write the centrifuge SOP for your clinic including balancing and cleaning.

✨ Grooming & professional tip

Photograph your prepared syringe next to the labelled tube for the record. It takes five seconds and settles any future question.

SELF-CHECK · సెల్ఫ్ చెక్

- ✦ Did I follow hygiene and consent for every client today?
- ✦ Which step did I need help with — and what will I drill tonight?
- ✦ Did I write my records and settings before leaving?

Trainer sign: _____ Trainee sign: _____